

# Taq DNA Polymerase (Mg<sup>2+</sup> free buffer)

Catalog # P102



Version 5.1

Vazyme biotech co., ltd.

## Introduction

Taq DNA Polymerase is a thermostable DNA polymerase that exhibits a 5'→3' polymerase activity and a 5'→3' exonuclease activity, with no 3'→5' exonuclease activity. Taq DNA Polymerase is purified from an *Escherichia coli* (*E. coli*) strain overexpressing the gene of *Thermus aquaticus* DNA Polymerase. No endonuclease, exonuclease, or bacterial DNA were detected in this kit.

The PCR products contain A at the 3'-end and can be directly cloned into T-Vectors. The products are also compatible with Vazyme ClonExpress II One Step Cloning Kit series (Vazyme, C#C112, #C113).

## Package Information

Components	P102-01 1,000 U	P102-d1 1,000 U	P102-02/d2 5,000 U	P102-03/d3 10,000 U
10x Taq Buffer (Mg <sup>2+</sup> free)	4 ml	4 ml		
25 mM MgCl <sub>2</sub>	4 ml	4 ml		
dNTP Mix (10 mM each)	-----	800 µl		
Taq DNA Polymerase (5 U/µl)	200 µl	200 µl	P102-01/d1 × 5	P102-01/d1 × 10

## Storage

Store at -20°C.

## Unit Definition

One unit (U) is defined as the amount of enzyme that incorporates 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C, with activated salmon sperm DNA as the template / primer.

## Quality Control

**Exonuclease Activity:** The product is tested in a reaction containing 10 U of Taq DNA Polymerase and 0.6 µg of λ-Hind III. After incubation at 37°C for 16 hours, there is no visually discernible change to DNA bands determined by agarose gel electrophoresis.

**Endonuclease Activity:** The product is tested in a reaction containing 10 U of Taq DNA Polymerase and 0.6 µg of Supercoiled pBR322. After incubation at 37°C for 4 hours, there is no visually discernible change to DNA band determined by agarose gel electrophoresis.

**Functional Assay:** The human α-1-antitrypsin gene is amplified for 30 cycles in a 50 µl system using 1.25 U of Taq DNA Polymerase and 100 ng human genomic DNA as template. A single DNA band of 360 bp is detected by 1% agarose gel electrophoresis.

## Protocol

### 1. General reaction mixture for PCR:

10x Taq Buffer (Mg <sup>2+</sup> free)	5 µl
25 mM MgCl <sub>2</sub> <sup>a</sup>	optional
dNTP Mix (10 mM each)	1 µl
Template DNA <sup>b</sup>	optional
Primer 1 (10 µM)	2 µl
Primer 2 (10 µM)	2 µl
Taq DNA Polymerase (5 U/µl) <sup>c</sup>	0.5 µl
ddH <sub>2</sub> O	to 50 µl

a. The final concentration of Mg<sup>2+</sup> of this mixture is 2 mM, as for most PCR reactions, the optimized volume of 25 mM MgCl<sub>2</sub> in a 50 µl PCR system is 3 µl-4 µl. However, if necessary, the concentration of Mg<sup>2+</sup> can be increased by adding 25 mM MgCl<sub>2</sub>.

b. The recommended amount of DNA template for a 50 µl reaction system is as follows:

Human Genomic DNA	0.1 - 1 µg
Bacterial Genomic DNA	10 - 100 ng
λ DNA	0.5 - 5 ng
Plasmid DNA	0.1 - 10 ng

c. The amount of Taq DNA Polymerase can be adjusted between 0.25 µl and 1 µl. Generally, higher level of Taq will increase the yield of PCR products but may decrease the specificity of PCR amplification.



Vazyme Biotech Co., Ltd  
www.vazyme.com

Order: global@vazyme.com

Support: support@vazyme.com

**For research use only, not for use in diagnostic procedures.**

## 2. Thermocycling conditions for a routine PCR:

---

94°C	5 min (Pre-denaturation)	} 30 - 35 cycles
94°C	30 sec	
55°C*	30 sec	
72°C	60 sec / kb	
72°C	7 min (Final extension)	

---

\*The optimal annealing temperature should be 1°C-2°C lower than the  $T_m$  of the primers used.

### Handling Notes

Taq DNA Polymerase also shows polymerase activity at room temperature. Therefore, it is recommended to set up reaction systems on ice and then immediately start the reaction in a PCR amplifier, so as to reduce nonspecific amplification during preparation and get better PCR results.

### Primers Designing Notes

1. Choose C or G as the last base of the 3'-end of the primer;
2. Avoid continuous mismatching at the last 8 bases of the 3'-end of the primer;
3. Avoid hairpin structure at the 3'-end of the primer;
4.  $T_m$  of the primers should be within the range of 55°C - 65°C;
5. Additional sequence should not be included when calculating  $T_m$  of the primers;
6. GC content of the primers should be within the range of 40% - 60%;
7.  $T_m$  and GC content of forward and reverse primers should be as similar as possible.

