# M-MLV (H-) Reverse Transcriptase



## Catalog # R021

Version 5.1

Vazyme biotech co., Itd.

#### Introduction

The wildtype Moloney Murine Leukemia Virus (M-MLV) reverse transcriptase has the following activities: RNA-dependent DNA polymerase, DNA-dependent DNA polymerase, and RNase H. The M-MLV (H-) Reverse Transcriptase is a single-site mutant of M-MLV which contains no RNase H activity. Compared with M-MLV mutants obtained via deletion of the RNase H domian, this product, which retains a complete protein structure and polymerase activities, can be used for the synthesis of longer cDNA or the preparation of cDNA library.

#### **Contents of Kits**

| Components                                  | R021 10,000 U |
|---|---------------|
| 5× RT Buffer                                | 500 μl        |
| M-MLV (H-) Reverse Transcriptase (200 U/μI) | 50 µI         |

#### **Storage**

All components should be stored at -20°C.

#### **Unit Definition**

One unit (U) is defined as the amount of enzyme that incorporates 1 nmol of dTTPs into acid-insoluble products in 10 min at 37°C with Poly(rA)-Oligo (dT) as the template / primer.

#### **Protocol**

1. Mix the following components in a RNase-free centrifuge tube:

| RNase free ddH₂O                            | to 20 µl   |
|---|--|
| 5× RT Buffer                                | 4 µl   |
| dNTP Mix (10 mM each)                       | 1 μΙ   |
| Oligo (dT) <sub>18</sub> (50 µM)            |  |
| or Random hexamers (50 ng/μl)               | 1 μΙ   |
| or Gene Specific Primers (2 μM)             |  |
| RNase inhibitor (40 U/µI)                   | 1 μΙ   |
| M-MLV (H-) Reverse Transcriptase (200 U/μl) | 1 μΙ   |
| Template RNA                                | Total RNA: 100 pg-5 µg Poly (A)⁺ RNA: 10 pg-500 ng |

### 2. Programs for the 1st-strand cDNA synthesis:

| For | oligo | (Th) |  |
|-----|-------|------|--|

70°C

| 42℃            | 45 min*      |  |
|----------------|--------------|--|
| 70℃            | 15 min       |  |
| For Random He  | examers      |  |
| 25℃            | 10 min       |  |
| 42℃            | 45 min*      |  |
| 70°C           | 15 min       |  |
| For Gene Speci | ific Primers |  |
| 42℃            | 45 min*      |  |

<sup>\*</sup> Can be optimized between 30 min and 60 min. Longer time is helpful to obtain longer cDNA (> 5 kb).

15 min

<sup>3.</sup> Incubate at 70°C for 15 min to inactivate the reverse transcriptase. The cDNA can be used for PCR or be stored at -20°C immediately. For PCR, it is recommended that the volume of cDNA  $\leq$  1/10 of total PCR reaction system volume.





Vazyme Biotech Co., Ltd Order: global@vazyme.com Support: support@vazyme.com For research use only, not for use in diagnostic procedures.