

T7 Endonuclease I

Catalog # EN303-01/02



Version 5.1

Vazyme biotech co., ltd.

1. Introduction

T7 Endonuclease I recognizes and cuts mismatched DNA, cruciform DNA, Holliday structures or junctions, and heteroduplex DNA. It can also recognize and cut double-strand DNA with nicks with a lower speed. The enzyme cuts the first, second or third phosphodiester bond at 5' end of the mismatch site. This product is a high purity T7 Endonuclease I protein with high activity, which stems from the expression of a recombinant T7 Endonuclease I (T7 Endonuclease I) gene in E. coli expression system.

2. Package Information

Components	EN303-01 (250U)	EN303-02 (1,250U)
T7 Endonuclease I	25 μ l	125 μ l
T7 Endonuclease I Reaction Buffer(10x)	1 ml	1 ml

3. Storage

Store at -20°C.

4. Reaction conditions

1 \times T7 Endonuclease I Reaction Buffer, incubate at 37°C.

5. Unit definition

1 unit is defined as the amount of enzyme required to cut more than 90% of 1 μ g of supercoiled cross-structured pUC (AT) plasmid into more than 90% of linearized DNA at 37°C within 1 hour in a 50 μ l reaction system.

6. Applications

Detection of gene mutation, SNP, mutations produced by TALEN and CRISPR/CAS9.
Recognition of mismatched DNA, resolve four-way junction or branched DNA.
Detecting or cutting heterogeneous dimer of DNA.
Random cutting of linear DNA for shot-gun cloning.

7. Strategy for detection of mutants with T7 Endonuclease I

1. PCR amplification of DNA fragments with mutation sites

1. Extract genomic DNA of cells after transfection.
2. DNA fragments with mutation sites, such as the target sites of TALEN or Cas9, are amplified through PCR, which should be performed with high fidelity DNA polymerase. The recommended size of the amplified fragment is 0.5-1 kb. Please do not set the mutation point in the middle of the amplified fragment, in order to produce two fragments with obviously different sizes after the cutting. Three PCR reactions should be set up for each sample with different templates as follows:
 - a: Genomic DNA of target cells (such as TALEN or Cas9 transfected cells)
 - b: Genomic DNA of negative control cells (cells without transfection of specific target sequence)
 - c: Water (the non-template control)
3. After the PCR reaction, do agarose gel electrophoresis analysis with a small amount of the PCR product. If the amplified products showed unique band with correct size, perform the next experiments.
4. Purify the PCR products with beads or PCR product purification kits.
5. Quantify the purified PCR products.



Vazyme Biotech Co., Ltd
www.vazyme.com

Order: global@vazyme.com

Support: global@vazyme.com

For research use only, not for use in diagnostic procedures.

2. Digestion Reaction of T7 Endonuclease I.

1. Prepare the reaction mixture with the purified PCR products in order as follows:

No.	1	2
Target Cell PCR products (200 ng)	x μ l	0
Negative control PCR products (200 ng)	0	x μ l
10 x T7 Endonuclease I Reaction Buffer	2 μ l	2 μ l
Nuclease-free Water	To 19 μ l	To 19 μ l

2. Do denaturation and annealing with a PCR machine using the following program:

95°C	5 min
95-85°C	-2°C /sec
85-25°C	-0.1°C/sec
Hold	4°C

3. Add T7 Endonuclease I to the annealed PCR products:

Annealed PCR product	19 μ l
T7 Endonuclease I	1 μ l

4. Incubate at 37°C for 15 min.

5. After the reaction, add 1.5 μ l of 0.25 M EDTA to terminate the reaction.

6. After termination, the products of digestion can be analyzed directly by 2% agarose gel electrophoresis.

Results

PCR product of target cell	+	-	+	-	
PCR product of negative control cell	-	+	-	+	
T7 Endonuclease I	+	+	-	-	DL2000




Figure 1. Digestion products analyzed by 2% agarose gel electrophoresis

8. Notes

1. T7 Endonuclease I is a substrate-structure selective enzyme. The enzyme acts on different DNA substrates with different activity. Therefore, the dosage of the enzyme and reaction time must be adjusted according to the types of substrates.

2. The nonspecific nuclease activity will increase when the reaction temperature exceeds 42°C. Keep the reaction temperature under 55°C, otherwise the nuclease activity will decrease.

